

EE51100033

Cod. 08.4268.99

Soluzione tampone pH 4.01 ±0.02/25°C (rosso) ml 250



Descrizione

Caratteristiche:

- etichetta personalizzata CDL
- Certificato riferibile NIST

Dati Tecnici

pH	4,01
Volume (ml)	250
Accuratezza (° C)	±0.02/25
Colore	Rosso

Varianti

Codice	Prodotto	pH	Volume (ml)	Colore
08.4268.99	EE51100033	4,01	250	Rosso
08.4269.99	EE51100043	7,00	250	Verde
08.4270.99	EE51100133	4,01	500	Rosso

08.4271.99	EE51100143	7,00	500	Verde
------------	------------	------	-----	-------