

ART. 55 SMALTATI

Cod. 03.1900.50



MORTAIO COMPLETAMENTE SMALTATO D.MM 65 ML 70
HALDENWANGER 55/0a

Descrizione

Mortai internamente smaltati, con becco.

Dati Tecnici

Dimensioni (Ø mm)	65
Capacità (ml)	70

Varianti

Codice	Prodotto	Dimensioni (Ø mm)	Capacità (ml)
03.1900.50	ART. 55 SMALTATI	65	70
03.1900.53	MORTAI Ø mm 80	80	100
03.1900.55	MORTAI Ø mm 90	90	170
03.1900.60	MORTAI Ø mm 115	115	330
03.1900.62	MORTAI Ø mm 130	130	400
03.1900.64	MORTAI Ø mm 150	150	700
03.1900.67	MORTAI Ø mm 180	180	1000

