

MORTAI Ø mm 180

Cod. 03.1900.67



MORTAIO COMPLETAMENTE SMALTATO D.MM 180 ML1000
HALDENWANGER 55/8

Descrizione

Mortai internamente smaltati

Dati Tecnici

Dimensioni (Ø mm)	180
Capacità (ml)	1000

Varianti

Codice	Prodotto	Dimensioni (Ø mm)	Capacità (ml)
03.1900.50	ART. 55 SMALTATI	65	70
03.1900.53	MORTAI Ø mm 80	80	100
03.1900.55	MORTAI Ø mm 90	90	170
03.1900.60	MORTAI Ø mm 115	115	330
03.1900.62	MORTAI Ø mm 130	130	400
03.1900.64	MORTAI Ø mm 150	150	700
03.1900.67	MORTAI Ø mm 180	180	1000

